



Mouse Anti-Human VEGF monoclonal antibody, clone BK1 (CABT-L422M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant Human VEGF Protein.
Target	Human VEGF
Immunogen	Recombinant Human VEGF protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	BK1
Purification	Protein A or G purified from hybridoma culture supernatant
Conjugate	Unconjugated
Applications	ELISA(Cap), LFIA(Cap), LETIA(Cap). We recommend the following for sandwich ELISA (Capture - Detection): CABT-L422M - CABT-L423M
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS with 0.05% sodium azide
Preservative	0.05% sodium azide

Storage Store at -20°C.

Ship Wet ice

BACKGROUND

Keywords VEGFA; vascular endothelial growth factor A; VPF; VEGF; MVCD1; vascular permeability factor;
