



Mouse Anti-Cyclopiazonic acid monoclonal antibody, clone CPA (CABT-L4231)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	ELISA IC 50 (ppb): 3
Immunogen	Cyclopiazonic acid-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	CPA
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Wet ice

BACKGROUND

Introduction	Cyclopiazonic acid (CPA), an indole-tetramic acid, is a secondary metabolite produced by several genera of Aspergillus and Penicillium growing on stored grain, including sunflower seeds. ¹ It has been found in human and animal foods and food sources, including milk, eggs, and poultry
Keywords	Cyclopiazonic acid;CPA

GENE INFORMATION

Official Symbol	Cyclopiazonic acid
-----------------	--------------------