



Mouse Anti-Sulfa synergist agent monoclonal antibody, clone SSA (CABT-L4234)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	ELISA IC 50 (ppb): 1 Cross-Reactivity: Carvedilol 2ppb; TMP 1ppb; DVD 1ppb; OMP 2ppb
Immunogen	Sulfa synergist agent-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	SSA
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

BACKGROUND

Introduction

Sulfonamide or sulphonamide is the basis of several groups of drugs. The original antibacterial sulfonamides (sometimes called sulfa drugs or sulpha drugs) are synthetic antimicrobial agents that contain the sulfonamide group. Some sulfonamides are also devoid of antibacterial activity.

Keywords

Sulfonamides;1;3-Difluoro-5-methylsulfonylbenzene;3;5-Difluoro-1-(methylsulfonyl)benzene;2-Chloro-1-fluoro-4-methylsulfonylbenzene

GENE INFORMATION

Official Symbol

Sulfonamides
