



Mouse Anti-Robenidine monoclonal antibody, clone RBD (CABT-L4239)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	ELISA IC 50 (PPB): 0.5
Immunogen	Robenidine-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	RBD
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Wet ice

BACKGROUND

Introduction Robenidine is a coccidiostat. Robenidine is an antibiotic used for the control of coccidiosis, a debilitating protozoal infection in poultry. Although there are alternative antibiotics available, robenidine is important in the management of antibiotic resistance as farmers rotate the use of robenidine with other antibiotics to try to preserve the effectiveness of these products in fighting infections.

Keywords Robenidine;RBD;Robenidinum;Robenidina

GENE INFORMATION

Official Symbol Robenidine
