



Mouse Anti-Ethopabate monoclonal antibody, clone EPB (CABT-L4241)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	ELISA IC 50 (PPB): 250
Immunogen	Ethopabate-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	EPB
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Wet ice

BACKGROUND

Introduction	Ethopabate is a coccidiostat used in poultry
Keywords	Ethopabate;

GENE INFORMATION

Official Symbol	Ethopabate
-----------------	------------