



Mouse Anti-Triazophos monoclonal antibody, clone TZP (CABT-L4244)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	ELISA IC 50 (PPB): 5
Immunogen	Triazophos-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	TZP
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Wet ice

BACKGROUND

Introduction

Triazophos is an organic thiophosphate and an organothiophosphate insecticide. It has a role as an EC 3.1.1.7 (acetylcholinesterase) inhibitor, an acaricide, an agrochemical and a nematicide.

Keywords

Triazophos;Triazofos;O,O-Diethyl O-(1-phenyl-1H-1,2,4-triazol-3-yl) phosphorothioate

GENE INFORMATION

Official Symbol

Triazophos