



Mouse Anti-Isoprocarb monoclonal antibody, clone IPC (CABT-L4245)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	ELISA IC 50 (PPB): 10
Immunogen	Isoprocarb-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	IPC
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Wet ice

BACKGROUND

Introduction	Isoprocarb is a carbamate ester.
Keywords	Isoprocarb;IUPAC

GENE INFORMATION

Official Symbol	Isoprocarb
-----------------	------------