



Mouse Anti-Human IL-6 monoclonal Antibody, clone C21 (CABT-L4269)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	C21
Purification	Chromatography on protein G Sepharose
Conjugate	Unconjugated
Applications	ELISA(Det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-L4270 - CABT-L4269
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS, pH 7.4, 0.09 % sodium azide (NaN ₃)
Preservative	0.09% sodium azide
Storage	Store at 2-8 °C
Ship	Wet ice

BACKGROUND

Introduction

Interleukin-6 (IL-6) is a proinflammatory cytokine secreted by T cells and macrophages to stimulate immune response to trauma, especially burns or other tissue damage leading to inflammation. IL-6 is especially important in the early stages of T-cell differentiation. In this phase, it reinforces the effect of IL-2 and promotes the differentiation of CD4 cells into T-helper2 cells.

Keywords

IL6;interleukin 6;HGF;HSF;BSF2;IL-6;IFNB2;interleukin-6;CDF;BSF-2;IFN-beta-2;interferon beta-2;interleukin BSF-2;interferon;beta 2;hybridoma growth factor;CTL differentiation factor;B-cell stimulatory factor 2;B-cell differentiation factor;PF-04236921