



Mouse Anti-Human Fetal Fibronectin monoclonal antibody, clone W4C2 (CABT- L4279)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	W4C2
Purification	Affinity purified
Conjugate	Unconjugated
Applications	ELISA(Det), LFIA(Det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-L4278 - CABT-L4279
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	0.2µM filtered, supplied liquid in Phosphate Buffer Saline buffer at pH 7.0 containing 0.09% sodium azide.
Preservative	0.09% sodium azide
Storage	Store at -20° C for long term storage. Avoid Freeze/Thaw Cycles.
Ship	Wet ice

BACKGROUND

Introduction	Fetal fibronectin (fFN) is a fibronectin protein produced by fetal cells. It is found at the interface of the chorion and the decidua (between the fetal sac and the uterine lining). It can be thought of as an adhesive or "biological glue" that binds the fetal sac to the uterine lining.
Keywords	fibronectin protein;fFN;Fetal Fibronectin;Fetal FN;FN;CIG;FNZ;MSF;ED-B;FINC;GFND;LETS;GFND2