



# Mouse Anti-Human PD-1 (CD279) Monoclonal antibody, clone J110 (CABT-L4466)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

### Product Overview

The J110 monoclonal antibody reacts with human PD-1 (programmed death-1) also known as CD279. PD-1 is a 50-55 kDa cell surface receptor encoded by the *Pdcd1* gene that belongs to the CD28 family of the Ig superfamily. PD-1 is transiently expressed on CD4 and CD8 thymocytes as well as activated T and B lymphocytes and myeloid cells. PD-1 expression declines after successful elimination of antigen. Additionally, *Pdcd1* mRNA is expressed in developing B lymphocytes during the pro-B-cell stage. PD-1's structure includes a ITIM (immunoreceptor tyrosine-based inhibitory motif) suggesting that PD-1 negatively regulates TCR signals. PD-1 signals via binding its two ligands, PD-L1 and PD-L2 both members of the B7 family. Upon ligand binding, PD-1 signaling inhibits T-cell activation, leading to reduced proliferation, cytokine production, and T cell death. Additionally, PD-1 is known to play key roles in peripheral tolerance and prevention of autoimmune disease in mice as PD-1 knockout animals show dilated cardiomyopathy, splenomegaly, and loss of peripheral tolerance. Induced PD-L1 expression is common in many tumors including squamous cell carcinoma, colon adenocarcinoma, and breast adenocarcinoma. PD-L1 overexpression results in increased resistance of tumor cells to CD8 T cell mediated lysis. In mouse models of melanoma, tumor growth can be transiently arrested via treatment with antibodies which block the interaction between PD-L1 and its receptor PD-1. For these reasons anti-PD-1 mediated immunotherapies are currently being explored as cancer treatments.

<b>Target</b>	Human PD-1 (CD279)
<b>Immunogen</b>	CBA mouse thymocytes
<b>Isotype</b>	IgG1, κ
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	J110

<b>Purification</b>	Protein G purified. Purity>95%. Determined by SDS-PAGE
<b>Conjugate</b>	Functional Grade
<b>Applications</b>	FC
<b>Molecular Weight</b>	150 kDa
<b>Format</b>	0.2 µM filtered liquid. Purified from tissue culture supernatant in an animal free facility
<b>Concentration</b>	Lot specific
<b>Size</b>	5 mg
<b>Buffer</b>	PBS, pH 7.0. Contains no stabilizers or preservatives. [low endotoxin azide-free]  Endotoxin level: <2EU/mg (<0.002EU/µg). Determined by LAL gel clotting assay Related dilution buffer: CABT-LB04
<b>Preservative</b>	None
<b>Storage</b>	The antibody solution should be stored undiluted at 4°C, and protected from prolonged exposure to light. Do not freeze.
<b>Ship</b>	Wet ice

## BACKGROUND

<b>Introduction</b>	This gene encodes a cell surface membrane protein of the immunoglobulin superfamily. This protein is expressed in pro-B-cells and is thought to play a role in their differentiation. In mice, expression of this gene is induced in the thymus when anti-CD3 antibodies are injected and large numbers of thymocytes undergo apoptosis. Mice deficient for this gene bred on a BALB/c background developed dilated cardiomyopathy and died from congestive heart failure. These studies suggest that this gene product may also be important in T cell function and contribute to the prevention of autoimmune diseases. [provided by RefSeq, Jul 2008]
<b>Keywords</b>	PD1;programmed cell death 1;PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1;programmed cell death protein 1;protein PD-1;systemic lupus erythematosus susceptibility 2;

## GENE INFORMATION

<b>Official Symbol</b>	programmed cell death 1
<b>Synonyms</b>	PD1; programmed cell death 1; PD1; PD-1; CD279; SLEB2; hPD-1; hPD-I; hSLE1;

