



Mouse Anti-Human PINP monoclonal antibody, clone 151 (CABT-L4571)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes native PINP
Immunogen	Native PINP
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	151
Purification	Protein A affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, CLIA We recommend the following for sandwich ELISA (Capture - Detection): CABT-L4570 - CABT-L4571 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Preservative	PBS with Proclin 300, pH7.4
Storage	Store -20°C for long term.

BACKGROUND

Introduction	Amino-terminal propeptide of type I procollagen (PINP) is released into blood circulation during bone formation. PINP is used as a bone turnover marker for the assessment of fracture risk and monitoring of osteoporosis treatment.
Keywords	P1NP;Procollagen Type I Intact N-Terminal Propeptide
