



Mouse Anti-Hyaluronic acid binding protein (HABP) monoclonal antibody, clone CYE112 (CABT-L4621)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	HABP
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	CYE112
Purification	Caprylic acid ammonium sulfate purified. Purity>90% by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA, CLIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Wet ice

BACKGROUND

Introduction	Hyaluronic Acid Binding Protein (HABP) is a reported synonym of the human protein 'hyaluronan binding protein 2', encoded by the gene HABP2.
Keywords	HABP; hyaluronan binding protein;