



Mouse Anti-Scl-70 monoclonal antibody, clone 28G3 (CABT-L4628)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Chromatically purified mouse monoclonal antibody to the Y-F-V-P-P-K-E-D-F-K epitope of the Scl-70 antigen
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	28G3
Purification	This material has been purified from mouse ascites fluid using chromatographic techniques
Conjugate	Unconjugated
Applications	For use in standard solid phase immunoassays.
Reconstitution	2 ml Deionized water.
Format	Lyophilized
Size	200 µg
Buffer	20 mM Sodium Phosphate, 0.02% NaN ₃ , pH 7.0
Preservative	0.02% NaN ₃
Storage	Store frozen at –70°C or lower
Ship	Wet ice

BACKGROUND

Keywords

Scl-70; TOP1
