



# Rabbit Anti-Human TRIM74 (aa 47-249) polyclonal antibody (CABT-L480R)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant Human Tripartite motif-containing protein 74 protein. aa 47-249
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	>95%, Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Suitable for use in ELISA and IHC Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Buffer</b>	50% Glycerol, 0.01M PBS, PH 7.4 with 0.03% Proclin 300
<b>Storage</b>	Store at -20°C.
<b>Ship</b>	Wet ice

## BACKGROUND

**Keywords**

TRIM74; tripartite motif containing 74; TRIM50C; tripartite motif-containing protein 74; tripartite motif-containing 74; tripartite motif-containing 50C; tripartite motif-containing protein 50C;

---