



Rabbit Anti-Human DERL3 polyclonal antibody (CABT-L528R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Human DERL3
Immunogen	Synthetic peptide conjugated to KLH. Corresponding to amino acids 51-150 of Human DERL3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, IHC, ICC/IF
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Preservative	0.03% Proclin300
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Keywords

DERL3; Der1-like domain family, member 3; C22orf14, chromosome 22 open reading frame 14; derlin-3; derlin 3; FLJ43842; IZP6; MGC71803; DERtrin 3; DERtrin-3; der1-like protein 3; degradation in endoplasmic reticulum protein 3; LLN2; C22orf14;
