



Rabbit Anti-Human ATRX polyclonal antibody (CABT-L533R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal to ATRX - C-terminal
Target	Human ATRX
Immunogen	This ATRX antibody was raised against a peptide within the C-terminus of human ATRX.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	Purified IgG in 70 mM Tris (pH 8), 105 mM NaCl, 31 mM glycine, 0.07 mM EDTA, 30% glycerol and 0.035% sodium azide. Sodium azide is highly toxic.
Preservative	0.035% sodium azide
Storage	Store at -20°C.

BACKGROUND

Keywords

ATRX; alpha thalassemia/mental retardation syndrome X-linked; alpha thalassemia/mental retardation syndrome X linked (RAD54 (S. cerevisiae) homolog) , JMS, Juberg Marsidi syndrome , RAD54; transcriptional regulator ATRX; RAD54 homolog (S. cerevisiae); XH
