



Mouse Anti-HCV core monoclonal antibody, clone IN652 (CABT-L541MH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	HCV core (33-38 aa epitope), reactive with genotypes 1a, 1b & 2a
Target	HCV core
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	HCV
Clone	IN652
Purification	Affinity purified, >90% purity
Conjugate	Unconjugated
Applications	ELISA(Det), LFIA, IF, WB, IHC We recommend the following for sandwich ELISA(Capture - Detection): CABT-L540MH - CABT-L541MH
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Ship

Wet ice
