



Rabbit Anti-COL17A1 monoclonal antibody, clone TS57-16 (CABT-L575)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Collagen XVII
Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	TS57-16
Purification	Protein A purified.
Conjugate	Unconjugated
Applications	WB, ICC/IF, IHC, IP
Molecular Weight	130/180 kDa
Cellular Localization	Cell junction, Membrane, Secreted
Positive Control	HUVEC, mouse skin tissue.
Format	Liquid
Size	100 µl
Buffer	1×TBS (pH7.4), 1% BSA, 40% Glycerol.
Preservative	0.05% Sodium Azide

Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
----------------	--

BACKGROUND

Introduction	The extensive collagen family comprises several chain types, including fibril-forming interstitial collagens and basement membrane collagens, with each type containing multiple isoforms. Products of the COL gene family, collagens are characterized as fibrous, extracellular matrix proteins with high tensile strength that constitute the major components of connective tissues, such as tendons and cartilage. All collagens contain a triple helix domain and frequently show lateral self-association in order to form complex connective tissues. Collagen Type XVII, also designated BP180, represents a type II transmembrane, epithelial adhesion molecule that plays a role in cell migration and differentiation. The full length Collagen Type XVII protein is expressed in hemidesmosomes of keratinocytes. Proteolytic shedding of Collagen Type XVII results in a species in the extracellular matrix, and this process may be mediated by a disintegrin and metalloprotease (ADAM) family member. The BPAG2 gene, which encodes the Collagen Type XVII protein, maps to human chromosome 10q25.1. Mutations in this gene result in Bullous pemphigoid, an inflammatory subepidermal blistering skin disease associated with an IgG autoimmune response to Collagen Type XVII.
---------------------	--

Keywords	180 kDa bullous pemphigoid antigen 2;Alpha 1 type XVII collagen;BA16H23.2;BP 180;BP180;BPA 2;BPAG 2;BPAG2;Bullous pemphigoid antigen 2;COL17A1;Collagen 17;Collagen alpha 1 XVII chain;Collagen alpha 1(XVII) chain;Collagen alpha1 XVII chain;Collagen type XVII alpha 1;Collagen XVII alpha 1 polypeptide;CollagenXVII;Epidermolysis bullosa junctional localisata variant;FLJ60881;KIAA0204;LAD 1;LAD1 antibody
-----------------	---