



Rabbit Anti-p. gingivalis RgpB (arginine-gingipain B) polyclonal Antibody (CABT-L587R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Porphyromonas gingivalis Gingipain R2 protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Porphyromonas gingivalis
Purification	Affinity purification
Conjugate	Unconjugated
Applications	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Buffer	50% Glycerol, 0.01M PBS, pH 7.4 with 0.03% Proclin 300
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Keywords

gpB; prtRII; rgp2; PG_0506Gingipain R2; EC 3.4.22.37; Arg-gingipain; Gingipain 2; RGP-2
