



# Rabbit Anti-Human CENPA polyclonal antibody (CABT-L599R)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CENPA
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IF Recommended dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Preservative</b>	0.02% sodium azide
<b>Storage</b>	Store at -20°C.
<b>Ship</b>	Wet ice

## BACKGROUND

**Keywords**

CENPA; centromere protein A; Cenp-A; histone H3-like centromeric protein A; centrosomin A; centromere autoantigen A;

---