



Rabbit Anti-GFP monoclonal antibody, clone TQ170-1 (CABT-L600)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	GFP
Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Clone	TQ170-1
Purification	Protein A purified.
Conjugate	Unconjugated
Applications	WB, ICC/IF, IHC
Positive Control	K562 cells transfected with GFP.
Format	Liquid
Size	100 µl
Buffer	1×TBS (pH7.4), 1% BSA, 40% Glycerol.
Preservative	0.05% Sodium Azide
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

The green fluorescent protein (GFP) was originally identified as a protein involved in the bioluminescence of the jellyfish *Aequorea victoria*. GFP cDNA produces a fluorescent product when expressed in prokaryotic cells, without the need for exogenous substrates or cofactors, making GFP a useful tool for monitoring gene expression and protein localization in vivo. Several GFP mutants have been developed, including EGFP, which fluoresce more intensely than the wildtype GFP and have shifted excitation maxima, making them useful for FACS and fluorescence microscopy as well as double-labeling applications. GFP is widely used in expression vectors as a fusion protein tag, allowing expression and monitoring of heterologous proteins fused to GFP.

Keywords GFP antibodyGreen fluorescent protein antibodyyfp antibody
