



Mouse Anti-Human CENPA monoclonal antibody, clone 6B8 (CABT-L600M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a portion of human CENP-A
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	6B8
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, ICC/IF Recommended dilution: Western Blot 1:1000 1 ug/ml of CENP-A Antibody was sufficient for detection of CENPA in 20 ug of U2OS cell lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary Antibody.
Format	Liquid
Concentration	Lot specific
Buffer	PBS, 50% Glycerol
Preservative	0.09% Sodium Azide
Storage	Store at -20°C.

Ship

Wet ice

BACKGROUND

Keywords

CENPA; centromere protein A; Cenp-A; histone H3-like centromeric protein A; centrosomin A; centromere autoantigen A;
