



# Rabbit Anti-Chicken RBC (Red Blood Cell) polyclonal antibody (Antiserum) (CABT-L6057)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Chicken RBC
Isotype	Serum
Source/Host	Rabbit
Species Reactivity	Chicken
Purification	Unpurified
Conjugate	Unconjugated
Applications	Hemagglutination Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Size	2 ml
Buffer	Antiserum
Preservative	0.05% sodium azide
Storage	Divide into useful aliquots and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

## BACKGROUND

**Keywords**

RBC; Erythrocytes; RBCs; RBC; Red Blood Cells; Red Blood Cell

---