



Rabbit Anti-Human Osteopontin/OPN (59-74) polyclonal antibody (CABT-L6059)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human Osteopontin (aa 59-74), human and mouse Osteopontin. There was no cross reactivity obtained with human Osteonectin and human Bone Sialoprotein.
Immunogen	Synthetic human Osteopontin (aa 59-74) KLH-conjugated(PSKSNESHDMDDMDD)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Unpurified Antiserum
Conjugate	Unconjugated
Applications	ELISA, IHC, WB Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Reconstitution	Resuspend in 20 ul / 100 ul aqua bidest.
Format	Lyophilized.
Size	20 µl; 100 µl
Preservative	None
Storage	2-8°C (lyophilized); -20 °C(dissolved) Repeated thawing and freezing must be avoided

Ship

Wet ice
