



Rabbit Anti-Human Osteopontin/OPN (164-179) polyclonal antibody (CABT-L6060)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Specificity | Human Osteopontin (aa 164-179) There was no cross reactivity obtained with human Osteopontin peptide (aa 59-74), human Osteonectin, and human Bone Sialoprotein |
| Immunogen | Synthetic human Osteopontin (aa 164-179) KLH-conjugated(QYPDATDEDITSHMES) |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Unpurified Antiserum |
| Conjugate | Unconjugated |
| Applications | ELISA, IHC Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. |
| Reconstitution | Resuspend in 20 ul / 100 ul aqua bidest. |
| Format | Lyophilized. |
| Size | 20 µl; 100 µl |
| Preservative | None |
| Storage | 2-8°C (lyophilized); -20 °C(dissolved) Repeated thawing and freezing must be avoided |

Ship

Wet ice
