



Mouse Anti-Pig CD206 monoclonal antibody, clone 233E3 (CABT-L6066)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Mannose Receptor/CD206 was used as immunogen.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Pig, Human, Bovine, Sheep
Clone	233E3.
Purification	Ion exchange chromatography
Conjugate	Unconjugated
Applications	WB, FC, ICC/IF, IHC, IHC-F, IHC-P, IP Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded
Format	Liquid, Purified
Concentration	Lot specific
Size	100 µg
Buffer	Tris-NaCl (pH 8.0)
Preservative	None
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles

Ship

Wet ice
