



Rabbit Anti-TACC3 monoclonal antibody, clone TZ21-17 (CABT-L608)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	TACC3
Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	TZ21-17
Purification	Protein A purified.
Conjugate	Unconjugated
Applications	WB, ICC/IF, IHC
Molecular Weight	140 kDa
Cellular Localization	Cytoplasm.
Positive Control	HT29, Hela, HepG2, human tonsil tissue, human spleen tissue, mouse testis tissue.
Format	Liquid
Size	100 µl
Buffer	1×TBS (pH7.4), 1% BSA, 40% Glycerol.
Preservative	0.05% Sodium Azide

Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
---------	--

BACKGROUND

Introduction	TACC1 (transforming acidic coiled coil gene 1) is one of three TACC family members, which are thought to be involved in breast tumorigenesis. TACC1 is located on 8p11 chromosomal region that is amplified in approximately 15% of all breast tumor samples. The short arm of chromosome 8 also contains FGFR1 whose expression is enhanced in most breast cancer tumors. TACC family members, TACC1, TACC2, and TACC3, map very closely to the corresponding FGFR1, FGFR2, FGFR3 genes on chromosomes 8, 10, and 4. Subsequently, since they are phylogenetically related, it is proposed that TACC and FGFR have similar roles in cell growth and differentiation. Also, TACC1 contains a conserved C-terminal region as in the Drosophila homolog, D-TACC. It has been shown that D-TACC is necessary for normal spindle function, and the mammalian TACC proteins appears to interact with centrosomes and microtubules in a similar manner.
--------------	---

Keywords	ERIC 1;ERIC-1;ERIC1;MGC117382;MGC133242;OTTHUMP00000113796;TACC3;TACC3_HUMAN;Transforming acidic coiled coil containing protein 3;Transforming acidic coiled-coil-containing protein 3 antibody
----------	---