



Human Anti-Human TG monoclonal antibody, clone D16 (CABT-L6107)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human Anti-Human TG monoclonal antibody, clone D16
Isotype	IgG
Source/Host	Human
Clone	D16
Purification	Affinity purified.
Conjugate	Unconjugated
Applications	<p>ELISA, Control</p> <p>Critical components in assay development are standards or positive controls. This humanized antibody can be used as an alternative to serum based calibrators and controls.</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Liquid, Purified
Size	1 mg
Buffer	PBS or Tris
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
Ship	Wet ice