



# Mouse Anti-Human IL-1B monoclonal antibody, clone D6 (CABT-L6110)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Recognizes human IL1 Beta
Immunogen	Recombinant Protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	D6
Purification	Affinity purified
Conjugate	Unconjugated
Applications	ELISA, CLIA We recommend the following for sandwich ELISA (Capture - Detection): CABT-L6109 - CABT-L6110 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded
Format	Liquid, Purified
Concentration	Lot specific
Size	1 mg
Buffer	PBS containing 10% glycerol and 0.02% sodium azide.

<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
<b>Ship</b>	Wet ice