



Mouse Anti Diclazuril monoclonal antibody, clone DIC (CABT-L6123)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	DIC
Conjugate	Unconjugated
Applications	ELISA, LFIA Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid, Purified
Concentration	Lot specific
Size	100ug, 500ug, 1 mg
Buffer	PBS containing 10% glycerol and 0.02% sodium azide.
Preservative	None
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
Ship	Wet ice