



# Human Anti-Human alpha Synuclein monoclonal antibody, clone A1 (CABT-L6126)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	Human alpha Synuclein
Immunogen	Recombinant Protein
Isotype	IgG
Source/Host	Human
Species Reactivity	Human
Clone	A1
Purification	Protein A purified. Purity $\geq 95\%$ as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	<p>ELISA</p> <p>We recommend the following for sandwich ELISA (Capture - Detection):</p> <p>CABT-L6126 - CABT-L6127</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Liquid, Purified
Concentration	Lot specific
Size	1 mg
Buffer	PBS or Tris-Gly

<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
<b>Ship</b>	Wet ice