



# Human Anti-human IL-10 monoclonal antibody, clone 97 (CABT-L6169)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Human Anti-human IL-10 monoclonal antibody for ELISA
<b>Target</b>	Human IL-10
<b>Immunogen</b>	Recombinant Protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Human
<b>Species Reactivity</b>	Human
<b>Clone</b>	97
<b>Purification</b>	Protein A purified. Purity $\geq 95\%$ as determined by SDS-PAGE
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	<p>ELISA</p> <p>We recommend the following for sandwich ELISA (Capture - Detection):</p> <p>CABT-L6168 - CABT-L6169</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
<b>Format</b>	Liquid, Purified
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg

<b>Buffer</b>	PBS or Tris-Gly
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
<b>Ship</b>	Wet ice