



Human Anti-human IL-12/IL-23 p40 monoclonal antibody, clone 3H2 (CABT-L6172)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human Anti-human IL-12/IL-23 p40 monoclonal antibody for ELISA
Target	Human IL-12/IL-23 p40
Immunogen	Recombinant Protein
Isotype	IgG
Source/Host	Human
Species Reactivity	Human
Clone	3H2
Purification	Protein A purified. Purity $\geq 95\%$ as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA We recommend the following for sandwich ELISA (Capture - Detection): CABT-L6172 - CABT-L6173 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid, Purified
Concentration	Lot specific
Size	1 mg

Buffer	PBS or Tris-Gly
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
Ship	Wet ice