



Human Anti-human IL-2RA monoclonal antibody, clone AB2 (CABT-L6192)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human Anti-human IL-2RA monoclonal antibody for ELISA
Target	Human IL-2RA
Immunogen	Recombinant Protein
Isotype	IgG
Source/Host	Human
Species Reactivity	Human
Clone	AB2
Purification	Protein A purified. Purity $\geq 95\%$ as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA We recommend the following for sandwich ELISA (Capture - Detection): CABT-L6192 - CABT-L6193 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid, Purified
Concentration	Lot specific
Size	1 mg

Buffer	PBS or Tris-Gly
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
Ship	Wet ice