



# Human Anti-human IL-8 monoclonal antibody, clone 1F (CABT-L6204)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |  |
|---------------------------|--|
| <b>Product Overview</b>   | Human Anti-human IL-8 monoclonal antibody for ELISA  |
| <b>Target</b>             | Human IL-8   |
| <b>Immunogen</b>          | Recombinant Protein  |
| <b>Isotype</b>            | IgG  |
| <b>Source/Host</b>        | Human  |
| <b>Species Reactivity</b> | Human  |
| <b>Clone</b>              | 1F   |
| <b>Purification</b>       | Protein A purified. Purity $\geq 95\%$ as determined by SDS-PAGE   |
| <b>Conjugate</b>          | Unconjugated   |
| <b>Applications</b>       | <p>ELISA</p> <p>We recommend the following for sandwich ELISA (Capture - Detection):</p> <p>CABT-L6204 - CABT-L6205</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p> |
| <b>Format</b>             | Liquid, Purified   |
| <b>Concentration</b>      | Lot specific   |
| <b>Size</b>               | 1 mg   |

|                |   |
|----------------|---|
| <b>Buffer</b>  | PBS or Tris-Gly   |
| <b>Storage</b> | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles |
| <b>Ship</b>    | Wet ice   |