



Human Anti-human IGF-1 monoclonal antibody, clone 6E (CABT-L6227)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human Anti-human IGF-1 monoclonal antibody for ELISA
Target	Human IGF-1
Immunogen	Recombinant Protein
Isotype	IgG
Source/Host	Human
Species Reactivity	Human
Clone	6E
Purification	Protein A purified. Purity $\geq 95\%$ as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA We recommend the following for sandwich ELISA (Capture - Detection): CABT-L6227 - CABT-L6228 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid, Purified
Concentration	Lot specific
Size	1 mg

Buffer	PBS or Tris-Gly
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
Ship	Wet ice