



Human Anti-human AFP-L3 monoclonal antibody, clone 3C (CABT-L6250)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Product Overview | Human Anti-human AFP-L3 monoclonal antibody for ELISA |
| Specificity | Human AFP-L3. No cross reaction with AFP. |
| Target | Human AFP-L3 |
| Immunogen | Recombinant Protein |
| Isotype | IgG |
| Source/Host | Human |
| Species Reactivity | Human |
| Clone | 3C |
| Purification | Protein A purified. Purity $\geq 95\%$ as determined by SDS-PAGE |
| Conjugate | Unconjugated |
| Applications | ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. |
| Format | Liquid, Purified |
| Concentration | Lot specific |
| Size | 1 mg |
| Buffer | PBS or Tris-Gly |

| | |
|----------------|---|
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles |
| Ship | Wet ice |