



Human Anti-human AFP-L3 monoclonal antibody, clone 4D (CABT-L6251)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human Anti-human AFP-L3 monoclonal antibody for ELISA
Specificity	Human AFP-L3. No cross reaction with AFP.
Target	Human AFP-L3
Immunogen	Recombinant Protein
Isotype	IgG
Source/Host	Human
Species Reactivity	Human
Clone	4D
Purification	Protein A purified. Purity $\geq 95\%$ as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid, Purified
Concentration	Lot specific
Size	1 mg
Buffer	PBS or Tris-Gly

Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
Ship	Wet ice