



Human Anti-human GP73 monoclonal antibody, clone 4B3 (CABT-L6258)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human Anti-human GP73 monoclonal antibody for ELISA
Target	Human GP73
Immunogen	Recombinant Protein
Isotype	IgG
Source/Host	Human
Species Reactivity	Human
Clone	4B3
Purification	Protein A purified. Purity $\geq 95\%$ as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	<p>ELISA</p> <p>We recommend the following for sandwich ELISA (Capture - Detection):</p> <p>CABT-L6257 - CABT-L6258</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Liquid, Purified
Concentration	Lot specific
Size	1 mg

Buffer	PBS or Tris-Gly
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
Ship	Wet ice