



Chicken Anti-VZV Polyclonal Antibody (CABT-L6331)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	VZV
Immunogen	Partially purified infected cell extract strain Ellen
Isotype	IgY
Source/Host	Chicken
Species Reactivity	VZV
Purification	Protein A purified
Conjugate	Unconjugated
Applications	IFA, IP, ELISA, WB Reactive with VZV in IFA, IP, ELISA, WB at 10 µg/ml. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	2 mg/ml
Size	1 mL
Buffer	PBS pH 7.4
Preservative	None
Storage	This product is supplied frozen on dry ice. Upon receipt, store at -20°C. Avoid multiple freeze-

thaw cycles as product degradation may result.

Ship

Wet ice
