



Rabbit Anti-Human SNCA (21-70) polyclonal antibody (CABT-L6360)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Species reactivity: Human, Mouse, Rat
Immunogen	Synthesized peptide derived from the N-terminal region of Human Synuclein- α .
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	Each laboratory should determine an optimum working titer for use in its particular application.
Positive Control	Synthesized peptide derived from the N-terminal region of Human Synuclein- α .
Format	Liquid
Concentration	Lot specific
Size	50 ug, 100 ug
Buffer	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
Preservative	0.02% sodium azide
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Ship	Wet ice

BACKGROUND

Keywords

SNCA; synuclein, alpha (non A4 component of amyloid precursor); PD1; NACP; PARK1;
PARK4; alpha-synuclein; synuclein alpha-140; non A-beta component of AD amyloid;
