



Rabbit Anti Human GALNT9 polyclonal antibody (CABT-L6364)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Human Polypeptide N-acetylgalactosaminyltransferase 9 protein (184-338AA)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	Recombinant Human Polypeptide N-acetylgalactosaminyltransferase 9 protein (184-338AA)
Format	Liquid
Concentration	Lot specific
Size	50 ug, 100 ug
Buffer	50% Glycerol, 0.01M PBS, pH 7.4, 0.03% Proclin 300
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Ship	Wet ice

BACKGROUND

Keywords GALNT9; UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylgalactosaminyltransferase 9 (GalNAc-T9); polypeptide N-acetylgalactosaminyltransferase 9; GALNAC T9; pp-GaNTase 9;

GalNAc transferase 9; polypeptide GalNAc transferase 9; protein-UDP
acetylgalactosaminyltransferase 9; UDP-GalNAc:polypeptide N-
acetylgalactosaminyltransferase 9; UDP-GalNAc: polypeptide N-
acetylgalactosaminyltransferase 9
