



# Rabbit Anti-Human VWF (aa 1-273) polyclonal antibody (CABT-L6476Z)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant Human von Willebrand factor protein (1-273AA)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC 1:20-1:200 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	50 µg, 100 µg
<b>Buffer</b>	50% Glycerol, 0.01M PBS, PH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

## GENE INFORMATION

UniProt ID

[P04275](#)

---