



Mouse Anti-Human β -NGF monoclonal antibody, clone N96 (CABT-L6521Z)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human β -NGF
Immunogen	E.coli derived Recombinant Human β -NGF
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human, Mouse
Clone	N96
Purification	Affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB. We recommend the following for sandwich ELISA (Capture - Detection): CABT-L6521Z - CABT-L6520Z
Preparation	Produced in mice immunized with highly pure Recombinant β -NGF. Anti- β -NGF-specific antibody was purified from cell culture by Protein G affinity chromatography
Format	Liquid
Concentration	Lot specific
Size	500 μ g
Buffer	PBS with preservation

Preservative	NaN3
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Ship	Wet ice