



Rabbit Anti-ITGAV monoclonal antibody, clone TD67-18 (CABT-L690)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Integrin alpha V
Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	TD67-18
Purification	Protein A purified.
Conjugate	Unconjugated
Applications	WB, FC
Molecular Weight	86/116 kDa
Cellular Localization	Membrane.
Positive Control	A549, PC-12.
Format	Liquid
Size	100 µl
Buffer	1×TBS (pH7.4), 1% BSA, 40% Glycerol.
Preservative	0.05% Sodium Azide

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Keywords

CD51 antibodyIntegrin alpha-V light chain antibodyintegrin, alpha V (vitronectin receptor, alpha polypeptide, antigen CD51) antibodyIntegrin, alpha V antibodyIntegrin, beta 6 antibodyITAV_HUMAN antibodyITGAV antibodyITGB6 antibodyVitronectin receptor subunit alpha antibody
