



Rabbit Anti-Human TSPYL5 Polyclonal antibody (CABT-L703R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide derived from human TSPYL5 (aa 251-350)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	Recommended dilution: ELISA: 1:500-1000; IHC: 1:100-500; IF: 1:100-500
Format	Liquid
Concentration	Lot specific
Size	200 µl
Buffer	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Preservative	None
Storage	Maintain at -20°C for up to 12 months. Avoid repeated freeze-thaw cycles. Store product undiluted.
Ship	Wet ice

BACKGROUND

Introduction

TSPYL5 (TSPY-like 5) is a protein-coding gene. Diseases associated with TSPYL5 include lung squamous cell carcinoma, and acute myeloid leukemia. GO annotations related to this gene include protein binding. An important paralog of this gene is SET

Keywords

TSPYL5; testis-specific protein, Y-encoded-like 5; 6330530B20; E130308C19Rik; testis-specific Y-encoded-like protein 5; TSPY-like protein 5;

GENE INFORMATION

Entrez Gene ID

[85453](#)

UniProt ID

[Q86VY4](#)