



Rabbit Anti-Human IRAKM monoclonal antibody, clone TO75-18 (CABT-L729)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	IRAKM
Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	TO75-18
Purification	Protein A purified.
Conjugate	Unconjugated
Applications	WB, IP
Molecular Weight	61 kDa
Cellular Localization	Cytoplasm, Nucleus.
Positive Control	Human placenta tissue, human thyroid tissue.
Format	Liquid
Size	100 µl
Buffer	1×TBS (pH7.4), 1% BSA, 40% Glycerol.
Preservative	0.05% Sodium Azide

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Interleukin-1 receptor (IL1R)-associated kinases (IRAKs) are important mediators in the signal transduction of Toll-like receptor (TLR) and IL1R family members, collectively referred to as TIRs. Binding of IL-1 to its cognate receptor results in the activation of the NFκB signaling pathway. An IL-1-dependent kinase termed IRAK-1 (for IL-1 receptor- associated kinase 1) coimmunoprecipitates with activated IL-1RI and is implicated as an upstream mediator of NFκB activation. A related Drosophila protein, Pelle, is a known upstream activator of Dorsal, the Drosophila homolog of NFκB. IRAK-2 is a proximal mediator of IL1, a component of the IL1R signaling complex, and is required for IL1R-induced NFκB activation. IRAK-4, like IRAK-1 and Pelle, has auto- and cross-phosphorylation kinase activity. IRAK-4 is strongly expressed in kidney and is also found in lung, testis, small intestine, breast, liver, and placenta. In contrast to the other IRAKs that are expressed in most cell types, IRAK-M is restricted to monocytic cells. IRAK-M mRNA transcripts are found predominantly in PBL and the monocytic cell lines U937 and THP-1.

Keywords

ASRT5;IL 1 receptor associated kinase M;IL-1 receptor-associated kinase M;IL1 receptor associated kinase M;Interleukin 1 receptor associated kinase 3;Interleukin 1 receptor associated kinase M;Interleukin-1 receptor-associated kinase 3;IRAK M;IRAK-3;IRAK-M;Irak3;IRAK3_HUMAN antibody

GENE INFORMATION

Entrez Gene ID

[11213](#)

UniProt ID

[Q9Y616](#)
