



Mouse Anti-PEDV S1 monoclonal antibody, clone 24N (CABT-L743M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Porcine epidemic diarrhea virus Spike protein S1 (PEDV-S1). This antibody reactive with PEDV strains, non-reactive with TGEV; not suitable for neutralization
Target	PEDV
Immunogen	PEDV strain "CV777"
Isotype	IgG
Source/Host	Mouse
Clone	24N
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide
Storage	Store at -20°C.
Ship	Wet ice