



Rabbit Anti-BCL11B monoclonal Antibody, clone 25B6 (CABT-L768R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Epitope mapping has revealed that the antibody recognizes an antigen within the 150 amino acids located at the extreme amino terminus of the protein, CTIP2-(1–150).
Target	BCL11B
Immunogen	Human fusion protein.9
Isotype	IgG, λ
Source/Host	Rabbit
Species Reactivity	Guinea Pig, Rat, Zebrafish, Human, Cynomolgus Monkey, Mouse
Clone	25B6
Purification	Protein A affinity purified
Conjugate	Unconjugated
Applications	CHIP, IHC-P, IP, WB, FC, IF, IHC-Fr
Format	Purified, Liquid
Concentration	Lot specific(generally 1 mg/ml)
Size	200 µg, 1mg
Buffer	PBS with 0.02% Proclin 300.
Preservative	0.02% Proclin 300

Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	This chimeric rabbit antibody was made using the variable domain sequences of the original Rat IgG2a format, for improved compatibility with existing reagents, assays and techniques.
---------------------	--

GENE INFORMATION

UniProt ID	Q9C0K0
-------------------	------------------------