



Mouse Anti-Isofenphos-methyl monoclonal antibody, clone IFM-1 (CABT-L7791)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|--|
| Immunogen | Isofenphos-methyl with carrier protein |
| Isotype | IgG |
| Source/Host | Mouse |
| Species Reactivity | N/A |
| Clone | IFM-1 |
| Purification | Affinity Purified |
| Conjugate | Unconjugated |
| Applications | ELISA, LFIA |
| Format | Liquid |
| Concentration | Lot specific |
| Size | 1 mg |
| Buffer | PBS |
| Preservative | None |
| Storage | Store at -20°C. |
| Ship | Wet ice |

BACKGROUND

Introduction C₁₄H₂₂NO₄PS; Isofenphos-methyl is a broad spectrum, organophosphorus insecticide, used to control soil pests in grain and oil crops such as click beetles, mole cricket and grub.

Keywords Isofenphos-methyl; IFM
